

Weqas
Unit 6, Parc Tŷ Glas
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750
Email: contact@weqas.com



EXTERNAL
QUALITY
ASSESSMENT



INTERNAL
QUALITY
CONTROL



REFERENCE
MEASUREMENT
SERVICES



EDUCATION &
TRAINING

Weqas

GLOBAL PROVIDER OF QUALITY
IN DIAGNOSTIC MEDICINE

Weqas Connect Lab EQA Case Studies

Gareth Davies / Sam Jones

Weqas

Summary Report

Analyte:

Overall Performance Category:

Running PI analyte –sample scores:

Standard Report

Analyte results table:

Scoring Table:

Running PI scores:

Bias Chart (Absolute):

Bias Chart (Relative):

Precision Scores:

Sample Histograms:

Method Summary Data:

Weqas
Unit 6, Parc Tŷ Glas
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750
Email: contact@weqas.com

Weqas

GLOBAL PROVIDER OF QUALITY
IN DIAGNOSTIC MEDICINE



EXTERNAL
QUALITY
ASSESSMENT



INTERNAL
QUALITY
CONTROL



REFERENCE
MEASUREMENT
SERVICES




EDUCATION &
TRAINING

Case 1 – Haematinics - Ferritin

Programme: **Haematinics** • Distribution Code: **HA1124**
 Distribution Start: 18-Nov-2024 • Distribution End: 02-Dec-2024 • Report Issued: n/a • Report Status: n/a
 Requested By: gareth@weqas.com

This Distribution	
Your % Poor PI	6%
Median All Participant % Poor PI	0%
97.5 th Centile % Poor PI	16%



PI Ranges	
<1	Good
1-2	Acceptable
>2	Poor

PI Scores						
Location	Location 1					
Instrument Name	Instrument 1	Instrument 2	Instrument 3	Instrument 4	Instrument 5	Instrument 6
Instrument Serial #	SN0001	SN0002	SN0003	SN0004	SN0005	SN0006
Ferritin	Acceptable	Good	N/A	N/A	N/A	Poor
B12	Good	Good	N/A	N/A	N/A	Acceptable
Folate	Good	Good	Good	N/A	N/A	Good
Iron	N/A	N/A	Good	Good	Good	N/A
TIBC	N/A	N/A	Good	Good	Acceptable	N/A
Transferrin	N/A	N/A	Warning	Good	Acceptable	N/A
Transferrin Saturation	N/A	N/A	Acceptable	Acceptable	Poor	N/A
Overall % poor PI	0%	0%	5%	0%	13%	17%
Overall % Non-return	0%	0%	0%	0%	0%	0%

Programme: Haematincs • Distribution Code: HA1124
 Distribution Start: 18-Nov-2024 • Distribution End: 02-Dec-2024 • Report Issued: • Report Status: n/a

Participant Code: WQ#### • Location: Location 1 • Analyser Name: Instrument 6 • Serial #: SN0006
 Date Samples received: n/a • Date of Analysis: n/a • Operator Details: n/a • Storage Conditions: n/a


Analyte: Ferritin • Method: Roche Elecsys • Kit: Elecsys Ferritin (07027273190)

		Ferritin µg/L		
Your reported result (µg/L)		Method: Roche Elecsys	Instrument Model: cobas e 801	Overall
Sample 1	1282	Mean	1549.11	1534.40
		SD	86.24	83.52
		Uncertainty	16.070	21.311
		n	45	24
Sample 2	998	Mean	1168.22	1150.33
		SD	62.24	54.21
		Uncertainty	11.597	13.831
		n	45	24
Sample 3	15	Mean	15.83	15.01
		SD	0.82	0.79
		Uncertainty	0.152	0.200
		n	45	24
Sample 4	52	Mean	51.79	51.03
		SD	1.91	1.78
		Uncertainty	0.355	0.455
		n	45	24


Scoring

Sample	1	2	3	4	Overall Performance
Target (Method Mean)	1549.11	1168.22	15.83	51.79	
Weqas TAE	216.33	163.38	3.2	8.2	
PI	-2.47	-2.08	-0.52	0.05	Poor


Running PI



Bias (Absolute)



Bias (Relative %)







Precision

	HA0624	HA0724	HA0824	HA0924	HA1024	This distribution: HA1124
Sy.x	10.64	0.55	2.97	18.83	0.48	17.09
IS	1	0	0	3	0	2

Sy.x is the average deviation from the best fit line and is an index of scatter.

IS Score	Interpretation
0 to 10	Good
11 to 150	Acceptable to Warning level
> 150	Unacceptable (incl. Curvilinear Data)



Legend: All Results (grey), My Method (purple), My Instrument (blue), My Result (orange)

Method Summary

HA1124	Sample 1		Sample 2		Sample 3		Sample 4	
Method Name	Mean	CV%	Mean	CV%	Mean	CV%	Mean	CV%
Beckman (n=4)	1032.14	7.4	785.99	11.5	8.22	13.4	29.36	5.0
Roche cobas c (n=2)	1568.00	0.9	1198.00	0.4	14.45	6.1	49.30	0.0
Roche Elecsys (n=45)	1549.11	5.6	1168.22	5.3	15.83	5.2	51.79	3.7
Siemens Advia (n=1)	997.00	0.0	804.00	0.0	9.00	0.0	31.00	0.0
Siemens Atellica (n=1)	1290.00	0.0	952.00	0.0	9.00	0.0	31.00	0.0

There are no Weqas or Participant supplied comments for Ferritin for this distribution.

Summary Report

Analyte: Ferritin

Overall Performance Category: Poor

Running PI analyte –sample scores: mostly good, some acceptable, more poor than median

Standard Report

Analyte results table: samples 1 and 2 significant negative bias to method, instrument and overall mean. Samples 3 and 4 seem to compare well

Scoring Table: samples 1 and 2 negative red PI, samples 3 and 4 green PIs, overall poor.

Running PI scores: generally good no other poor in last 6 distributions

Bias Chart (Absolute): current results at higher concs low compared to method mean (target value, x axis) and outside instrument SDs, results in red area. previous samples at these concs good. Linear regression shows $y = 0.83x + 7.69$. -17% negative proportional bias, +7.69 positive constant bias.

Bias Chart (Relative): samples at higher concs in red, previous samples at these concs good

Precision Scores: $Sy.x$ 17, IS 2, IS 0-10 = good

Sample Histograms: samples 1 and 2 negative bias compared to others, samples 3 and 4 central to normal distribution

Method Summary Data: Roche Elecsys dominant method group

Weqas
Unit 6, Parc Tŷ Glas
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750
Email: contact@weqas.com

Weqas

GLOBAL PROVIDER OF QUALITY
IN DIAGNOSTIC MEDICINE



EXTERNAL
QUALITY
ASSESSMENT



INTERNAL
QUALITY
CONTROL



REFERENCE
MEASUREMENT
SERVICES




EDUCATION &
TRAINING

Case 2 – Lipids - Triglyceride

Programme: **Lipids** • Distribution Code: **LP1124**
 Distribution Start: 18-Nov-2024 • Distribution End: 02-Dec-2024 • Report Issued: n/a • Report Status: n/a
 Requested By: gareth@weqas.com

This Distribution	
Your % Poor PI	38%
Median All Participant % Poor PI	0%
97.5 th Centile % Poor PI	33%



PI Ranges	
<1	Good
1-2	Acceptable
>2	Poor

PI Scores		
Location	Location 1	
Instrument Name	Instrument 1	Instrument 2
Instrument Serial #	SN0001	SN0002
Cholesterol	Poor	Poor
Triglyceride	Poor	Poor
HDL Cholesterol	Acceptable	Good
LDL Cholesterol	Warning	Warning
Overall % poor PI	38%	38%
Overall % Non-return	0%	0%

Participant Code: WQ#### • Location: Location 1 • Analyser Name: Instrument 2 • Serial #: SN0002
 Date Samples received: 02-Dec-2024 • Date of Analysis: 02-Dec-2024 • Operator Details: n/a • Storage Conditions: -70 °C



Analyte: Triglyceride • Method: Roche TRIGL • Kit: TRIGL (20767107 322)

Precision

	LP0624	LP0724	LP0824	LP0924	LP1024	This distribution: LP1124
Sy.x	0.04	0	0.02	0.02		0.21
IS	29	0	1	0		35

Sy.x is the average deviation from the best fit line and is an index of scatter.

IS Score	Interpretation
0 to 10	Good
11 to 150	Acceptable to Warning level
> 150	Unacceptable (incl. Curvilinear Data)



Method Summary

LP1124	Sample 1	Sample 2	Sample 3	Sample 4				
Method Name	Mean	CV%	Mean	CV%	Mean	CV%	Mean	CV%
GPO (n=22)	2.731	2.5	4.576	2.7	2.773	2.6	0.841	6.3
NMR Spectroscopy (n=)	n/a	n/a	n/a	n/a	n/a	n/a	n/a	n/a
Roche TRIGL (n=66)	2.752	2.1	4.599	2.2	2.773	2.3	0.879	3.0
Vitros (n=3)	3.390	0.6	5.093	0.3	3.683	0.9	0.983	3.3


There are no Weqas or Participant supplied comments for Triglyceride for this distribution.

Triglyceride mmol/L						
	Your reported result (mmol/L)		Method: Roche TRIGL	Instrument Model: cobas c 501	Overall	Scoring Reference Value
Sample 1	3.25	Mean	2.752	2.783	2.751	2.64
		SD	0.058	0.081	0.084	
		Uncertainty	0.0089	0.0282	0.0084	
		n	66	15	91	
Sample 2	5.79	Mean	4.599	4.656	4.600	4.48
		SD	0.099	0.176	0.113	
		Uncertainty	0.0152	0.0569	0.0148	
		n	66	15	91	
Sample 3	3.64	Mean	2.773	2.794	2.777	2.65
		SD	0.064	0.094	0.070	
		Uncertainty	0.0098	0.0302	0.0091	
		n	66	15	91	
Sample 4	0.87	Mean	0.879	0.888	0.871	0.81
		SD	0.027	0.038	0.042	
		Uncertainty	0.0041	0.0122	0.0055	
		n	66	15	91	


Scoring

Sample	1	2	3	4	Overall Performance
Target (Reference Value)	2.64	4.48	2.65	0.81	
Weqas TAE	0.401	0.638	0.402	0.164	
PI	3.05	4.11	4.93	0.73	Poor


Running PI



Bias (Absolute)



Bias (Relative %)



Summary Report

Analyte: Triglyceride

Overall Performance Category: Poor

Running PI analyte –sample scores: no issues before this distribution but non return for LP1024

Standard Report

Analyte results table: samples 1, 2, 3 positive bias to method, instrument, overall mean and ref value.

Scoring Table: samples 1, 2, 3 positive red PI, overall performance poor

Running PI scores: no issues before this distribution but non return for LP1024. 3 red PIs this dist, 2 > 4

Bias Chart (Absolute): 3 samples significant positive bias, linear regression shows $y = 1.34x - 0.16$, +34% positive proportional bias.

Bias Chart (Relative): 3 samples significant positive bias

Precision Scores: $Sy.x$ 0.21, IS 35 (11-150 acceptable to warning), precision acceptable

Sample Histograms: samples 1, 2, 3 positive bias, sample 2 off the scale

Method Summary Data: Roche slightly positive to GPO method, scored ref values available

Weqas
Unit 6, Parc Tŷ Glas
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750
Email: contact@weqas.com

Weqas

GLOBAL PROVIDER OF QUALITY
IN DIAGNOSTIC MEDICINE



EXTERNAL
QUALITY
ASSESSMENT



INTERNAL
QUALITY
CONTROL




REFERENCE
MEASUREMENT
SERVICES



EDUCATION &
TRAINING

Case 3 – Urine Chemistry - Urate



PI Scores			
Location	Location 1		
Instrument Name	Instrument 1	Instrument 2	Instrument 3
Instrument Serial #	SN0001	SN0002	SN0003
Sodium	N/A	N/A	Good
Potassium	N/A	N/A	Good
Chloride	N/A	N/A	Acceptable
Urea	N/A	N/A	Acceptable
Creatinine	N/A	Good	Good
Calcium	N/A	N/A	Good
Phosphate	N/A	N/A	Good
Protein	N/A	Acceptable	Acceptable
Protein / Creatinine Ratio	N/A	Good	Good
Albumin	N/A	N/A	Good
Albumin / Creatinine Ratio	N/A	N/A	Good
Magnesium	N/A	N/A	Poor
Urate	N/A	N/A	Poor
Osmolality	Good	N/A	N/A
Overall % poor PI	0%	0%	10%
Overall % Non-return	0%	0%	0%


Analyte: Urate • Method: Uricase / Peroxidase • Kit: UA2 (08058750190)

Urate mmol/L					
	Your reported result (mmol/L)		Method: Uricase / Peroxidase	Instrument Model: cobas c 503	Overall
Sample 1	0.83	Mean	1.2434	1.0760	1.2434
		SD	0.0379	0.2416	0.0379
		Uncertainty	0.00768	0.17435	0.00768
		n	38	3	38
Sample 2	1.55	Mean	2.2941	2.0340	2.2941
		SD	0.0561	0.4779	0.0561
		Uncertainty	0.01138	0.34492	0.01138
		n	38	3	38
Sample 3	1.44	Mean	2.0485	1.7927	2.0485
		SD	0.0628	0.3466	0.0628
		Uncertainty	0.01274	0.25011	0.01274
		n	38	3	38


Scoring

Sample	1	2	3	Overall Performance
Target (Method Mean)	1.2434	2.2941	2.0485	
Weqas TAE	0.2264	0.3239	0.2979	
PI	-3.65	-4.59	-4.09	Poor


Running PI



Bias (Absolute)



Bias (Relative %)




Precision

	UC0624	UC0724	UC0824	UC0924	UC1024	This distribution: UC1124
Sy.x	0.02	0.01		0.08	0.03	0.05
IS	1	9		118	428	34

Sy.x is the average deviation from the best fit line and is an Index of scatter.

IS Score	Interpretation
0 to 10	Good
11 to 150	Acceptable to Warning level
> 150	Unacceptable (incl. Curvilinear Data)



Method Summary

UC1124	Sample 1		Sample 2		Sample 3	
Method Name	Mean	CV%	Mean	CV%	Mean	CV%
Uricase / Peroxidase (n=38)	1.2434	3.0	2.2941	2.4	2.0485	3.1

There are no Weqas or Participant supplied comments for Urate for this distribution.

Summary Report

Analyte: Urine urate

Overall Performance Category: poor

Running PI analyte: sample scores: some poor on most distributions

Standard Report

Analyte results table: significant negative results for all 3 samples compared to method mean, instrument mean lower than method mean but n=3

Scoring Table: all 3 samples negative red PI, overall performance poor

Running PI scores: same pattern for UC0924 and UC1124, otherwise all good

Bias Chart (Absolute): negative bias for all 3 samples, 3 previous samples similar pattern, others compare well to target values, instrument SDs wide but n only 3 so these results affecting this SD. Linear regression shows $y = 0.70x - 0.04$. 30% negative proportional bias

Bias Chart (Relative): all 3 samples significant negative bias, all in the red

Precision Scores: $Sy.x$ 0.05, IS 34 (11-150 acceptable to warning), precision acceptable

Sample Histograms: all 3 samples significant negative bias, do not compare to other participants

Method Summary Data: all participants in same method group