

Weqas  
Unit 6, Parc Tŷ Glas  
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750  
Email: [contact@weqas.com](mailto:contact@weqas.com)



GLOBAL PROVIDER OF QUALITY  
IN DIAGNOSTIC MEDICINE



EXTERNAL  
QUALITY  
ASSESSMENT



INTERNAL  
QUALITY  
CONTROL



REFERENCE  
MEASUREMENT  
SERVICES



EDUCATION &  
TRAINING

## EQA Report Interpretation

Gareth Davies / Sam Jones

# Overview

- Why multiple samples are important
- The Weqas Report - Statistical Indices
  - Imprecision
  - Inaccuracy
- How do we set Analytical Performance Specifications (APS) and Target Values
- Weqas Connect Reports
- Problem Solving and Cases
- Case Study Example

# Why multiple samples are important


Identifies components of both Imprecision and Inaccuracy

Identifies systematic errors

Assesses method linearity – required for ISO 15189

Powerful error detection tool

# Gaussian Distribution



# The Weqas Report

## Statistical Indices

- Imprecision:  $Sy.x$ ,  $r$ ,  $IS$
- Inaccuracy: Bias plot,  $y=mx+c$

# Interpretation of Imprecision

**Sy.x** is the deviation around the best fit line and is an indication of scatter. It is used as an index of within run imprecision. This is given in the units of the analyte in question.

Correlation coefficient (**r value**) is also an index of within run imprecision. The closer to 1.0 the value is, the better the precision (less scatter about the best fit line).

**IS score** is derived from the correlation coefficient ( $IS = (1-r)*10,000$ ).

'r' value	Imprecision score	
0.9990 to 1.0000	0 to 10	- Good
0.9850 to 0.9989	11 to 150	- Acceptable to Warning level
< 0.9850	> 150	- Unacceptable (including Curvilinear Data)

# Imprecision indices

Standard deviation

$$SD = \sqrt{\frac{\sum (\bar{y} - y)^2}{d.f.}}$$

where

$y$  = observed value

$\bar{y}$  = expected value

d.f. = degrees of freedom

Standard dev of residuals

$$S_{y.x} = \sqrt{\frac{\sum (\tilde{y} - y)^2}{d.f.}}$$


where

$\tilde{y}$  = the value on the  
line of best fit

# Imprecision — Scatter around the “best fit line”, $Sy.x$

Sodium (mmol/l)

Measurement of scatter



$r = 1.000$   
 $Sy.x = 0$

All reported results are on the “best fit” line – no scatter

$r = 0.9812$   
 $Sy.x = 4.1$

Reported results are scattered around the “best fit” line – standard deviation of 4.1 mmol/L



## Imprecision – how to calculate $S_{y.x}$


Target value	Lab result	Line of best fit	Deviation	
$x$	$y$	$\tilde{y}$	$y - \tilde{y} = dy.x$	$dy.x^2$
111	108	109.9	-1.9	3.62
123.5	128	123.3	4.7	21.97
135.7	136	136.4	-0.4	0.16
148	144	149.6	-5.6	31.30
160.3	166	162.8	3.2	10.31
<i>slope</i>	1.07		$\sum dy.x$	$\sum dy.x^2$
<i>int</i>	-9.17		0.0	67.36
			$d.f = n-1 = 4$	
	<b><math>S_{y.x} =</math></b>	<b><math>\sqrt{\sum dy.x^2 / d.f}</math></b>	<b><math>\sqrt{67.4/4}</math></b>	<b><math>4.10</math></b>

# Interpretation of Inaccuracy

Represents the point at which the line crosses the y axis.  
This is shown in the units of measurement of each analyte.  
If deemed significant, this indicates a constant error.


$$y = mx + c$$

Represents the slope. If deemed significant -  
indicates a proportional error.




## ERRORS IN ACCURACY

### 2. SYSTEMATIC PROPORTIONAL (POSITIVE BIAS)




### ERRORS IN ACCURACY

#### 3. SYSTEMATIC PROPORTIONAL (NEGATIVE BIAS)





### ERRORS IN ACCURACY 4. CONSTANT (POSITIVE BIAS)




## ERRORS IN ACCURACY


### 5. CONSTANT (NEGATIVE BIAS)











	Imprecision		Inaccuracy		
	Random	Curvilinear	Prop	Mixed	Constant
Slope	No	Yes/No	Yes	Yes	No
Intercept	No	Yes/No	No	Yes	Yes
$S_{y,x}$	Yes	Yes	No	No	No
$r$	Yes	Yes	No	No	No

Weqas  
Unit 6, Parc Tŷ Glas  
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750  
Email: [contact@weqas.com](mailto:contact@weqas.com)

# Weqas

GLOBAL PROVIDER OF QUALITY  
IN DIAGNOSTIC MEDICINE



EXTERNAL  
QUALITY  
ASSESSMENT



INTERNAL  
QUALITY  
CONTROL



REFERENCE  
MEASUREMENT  
SERVICES



EDUCATION &  
TRAINING

## How do we set Analytical Performance Specifications (APS) and Target Values?

# The Weqas Report

## Target values used in Statistical Analysis

Hierarchy




Reference values – used for bias plot & SDI / PI calculation

Method mean – used for SDI / PI calc if no ref value and  $n \geq 8$

Overall mean – used for SDI / PI calc and bias plot if no ref value and  $n < 8$


Instrument mean – on report for information only

# Hierarchy for setting APS




# Interlaboratory variation (state of the art)

## FT4



Relationship of analyte conc. to SD follows a polynomial equation for most analytes  
and is non linear therefore using fixed % TE is not appropriate

# “State of the art” v Biology



Biological goals not achievable. Although Performance has improved over time.



Weqas  
Unit 6, Parc Tŷ Glas  
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750  
Email: [contact@weqas.com](mailto:contact@weqas.com)



GLOBAL PROVIDER OF QUALITY  
IN DIAGNOSTIC MEDICINE



EXTERNAL  
QUALITY  
ASSESSMENT



INTERNAL  
QUALITY  
CONTROL



REFERENCE  
MEASUREMENT  
SERVICES



EDUCATION &  
TRAINING

## Weqas Connect Reports

# Standard Quantitative Report

Participant Code: WQ00000

#####  
#####  
#####  
#####  
#####

Weqas


Unit 8, Parc Ty Glas  
Llanishen  
Cardiff  
United Kingdom  
CF14 5DU  
Tel: 02920 314750  
Email: office@weqas.com  
Programme Organiser: Mrs Annette Thomas

Weqas

Programme: ED Toxicology • Distribution Code: SP1223  
Distribution Start: 04-Dec-2023 • Distribution End: 18-Dec-2023 • Report Issued: 27-Dec-2023 • Report Status: Final  
Requested By: #####

This Distribution	
Your % Poor PI	9%
Median All Participant % Poor PI	0%
97.5 <sup>th</sup> Centile % Poor PI	40%

Running PI analyte-sample Scores



PI Ranges	
<1	Good
1-2	Acceptable
>2	Poor

PI Scores			
Location	Hospital 1	Hospital 2	
Instrument Name	#####	#####	#####
Instrument Serial #	SN0001	SN0002	SN0003
Salicylate	Good	Warning	Poor
Paracetamol	Good	Good	Warning
Ethanol	Good	N/A	Good
Overall % poor PI	0%	13%	17%
Overall % Non-return	0%	0%	0%

PI Code	Meaning
N/A	Not enrolled for this analyte
NRR	Analyte enrolled but no results returned
N/S	This quantitative analyte not scored for this distribution

## \*\*\*\* Please note\*\*\*\*

The NQAAP expects 100% compliance on EQA returns. If you are enrolled for an analyte and have not returned a result (denoted by NRR in your PI table), this will be treated as a poor performance score. Please let us know as soon as possible if you are no longer providing this analyte as part of your diagnostic service.

## Distribution Specific Comment

There is no specific comment for this distribution.

No participant comments have been submitted for any instruments for this distribution.

Participant Code: WQ00000

#####  
#####  
#####  
#####  
#####



Unit 6, Parc Ty Glas  
Llanishen  
Cardiff  
United Kingdom  
CF14 5DU  
Tel: 02920 314750  
Email: office@weqas.com  
Programme Organiser: Mrs Annette Thomas

Programme: **ED Toxicology** • Distribution Code: **SP1223**  
Distribution Start: 04-Dec-2023 • Distribution End: 18-Dec-2023 • Report Issued: 27-Dec-2023 • Report Status: Final  
Requested By: #####

This Distribution

Your % Poor PI	9%
Median All Participant % Poor PI	0%
97.5 <sup>th</sup> Centile % Poor PI	40%

Running PI analyte-sample Scores

Sample	Good (%)	Acceptable (%)	Poor (%)	Non-returns (%)	Not Scored (%)
SP0723	85	15	0	0	0
SP0823	85	15	0	0	0
SP0923	70	30	0	0	0
SP1023	95	5	0	0	0
SP1123	90	10	0	0	0
SP1223	85	10	5	0	0
SP1223 Median	85	15	0	0	0

PI Ranges

<1	Good
1-2	Acceptable
>2	Poor

PI Scores			
Location	Hospital 1	Hospital 2	
Instrument Name	#####	#####	#####
Instrument Serial #	SN0001	SN0002	SN0003
Salicylate	Good	Warning	Poor
Paracetamol	Good	Good	Warning
Ethanol	Good	N/A	Good
Overall % poor PI	0%	13%	17%
Overall % Non-return	0%	0%	0%

PI Code	Meaning
N/A	Not enrolled for this analyte
NRR	Analyte enrolled but no results returned
N/S	This quantitative analyte not scored for this distribution

\*\*\*\*\* **Please note**\*\*\*\*\*

The NQAAP expects 100% compliance on EQA returns. If you are enrolled for an analyte and have not returned a result (denoted by NRR in your PI table), this will be treated as a poor performance score. Please let us know as soon as possible if you are no longer providing this analyte as part of your diagnostic service.

Distribution Specific Comment

There is no specific comment for this distribution.

*No participant comments have been submitted for any instruments for this distribution.*

Participant Code: WQ00000 - Location: Hospital 1 - Analyser Name: ##### - Serial #: SN0001  
Date Samples received: 05-Dec-2023 - Date of Analysis: 06-Dec-2023 - Operator Details: n/a - Storage Conditions: 2 - 8 °C


Analyte: Salicylate - Method: Roche Enzymatic - Kit: SALI

Salicylate mg/L					
	Your reported result (mg/L)		Method: Roche Enzymatic	Instrument Model: cobas c 501	Overall
Sample 1	160	Mean	160.1	161.1	159.9
		SD	4.5	4.0	6.1
		Uncertainty	0.93	1.67	1.03
		n	37	9	54
Sample 2	250	Mean	252.5	253.3	252.9
		SD	6.2	4.0	7.3
		Uncertainty	1.28	1.66	1.24
		n	37	9	54
Sample 3	705	Mean	704.6	703.9	713.8
		SD	23.1	18.9	32.1
		Uncertainty	5.18	8.33	5.78
		n	31	8	48
Sample 4	531	Mean	537.1	534.6	536.4
		SD	12.8	13.4	18.4
		Uncertainty	2.63	5.58	3.12
		n	37	9	54


### Scoring

Sample	1	2	3	4	Overall Performance
Target: (Method Mean)	160.1	252.5	704.6	537.1	
Weqas TAE	13	18.2	68.6	45.1	
PI	-0.02	-0.28	0.01	-0.27	Good


### Running PI



### Bias (Absolute)



### Bias (Relative %)




### Precision

	SP0723	SP0823	SP0923	SP1023	SP1123	This distribution: SP1223
Syx	8.28	6.26	2.05	1.72	22.32	3.62
IS	4	2	0	0	27	1

Syx is the average deviation from the best fit line and is an index of scatter.

IS Score	Interpretation
0 to 10	Good
11 to 150	Acceptable to Warning level
> 150	Unacceptable (Incl. Curvilinear Data)



### Method Summary

SP1223	Sample 1		Sample 2		Sample 3		Sample 4	
Method Name	Mean	CV%	Mean	CV%	Mean	CV%	Mean	CV%
Abbott Enzymatic (n=2)	159.8	1.4	258.0	0.7	760.7	1.0	557.4	1.2
CLS Enzymatic (n=1)	162.6	0.0	262.6	0.0	747.3	0.0	549.7	0.0
Emit tox Salicylic Acid (n=9)	162.8	6.0	256.3	4.3	744.2	5.1	539.4	6.0
Roche Enzymatic (n=37)	160.1	2.8	252.5	2.5	704.6	3.3	537.1	2.4
Siemens Enzymatic (n=5)	150.6	4.3	244.3	4.1	699.7	1.6	513.1	1.9

There are no Weqas or Participant supplied comments for Salicylate for this distribution.

**Programme:** ED Toxicology • **Distribution Code:** SP1223  
 Distribution Start: 04-Dec-2023 • Distribution End: 18-Dec-2023 • Report Issued: 27-Dec-2023 • Report Status: Final

**Participant Code:** WQ00000 • **Location:** Hospital 1 • **Analysers Name:** ##### • **Serial #:** SN0001  
 Date Samples received: 05-Dec-2023 • Date of Analysis: 06-Dec-2023 • Operator Details: n/a • Storage Conditions: 2 - 8 °C


**Analyte:** Salicylate • **Method:** Roche Enzymatic • **Kit:** SALI

Salicylate mg/L					
	Your reported result (mg/L)		Method: Roche Enzymatic	Instrument Model: cobas c 501	Overall
Sample 1	160	Mean	160.1	161.1	159.9
		SD	4.5	4.0	6.1
		Uncertainty	0.93	1.67	1.03
		n	37	9	54
Sample 2	250	Mean	252.5	253.3	252.9
		SD	6.2	4.0	7.3
		Uncertainty	1.28	1.66	1.24
		n	37	9	54
Sample 3	705	Mean	704.6	703.9	713.8
		SD	23.1	18.9	32.1
		Uncertainty	5.18	8.33	5.78
		n	31	8	48
Sample 4	531	Mean	537.1	534.6	536.4
		SD	12.8	13.4	18.4
		Uncertainty	2.63	5.58	3.12
		n	37	9	54


## Scoring

Sample	1	2	3	4	Overall Performance
Target: (Method Mean)	160.1	252.5	704.6	537.1	
Weqas TAE	13	18.2	68.6	45.1	
PI	-0.02	-0.28	0.01	-0.27	Good


## Running PI



### Bias (Absolute)



### Bias (Relative %)









## Precision

	SP0723	SP0823	SP0923	SP1023	SP1123	This distribution: SP1223
Sy.x	8.28	6.26	2.05	1.72	22.32	3.62
IS	4	2	0	0	27	1

*Sy.x is the average deviation from the best fit line and is an index of scatter.*

IS Score	Interpretation
0 to 10	Good
11 to 150	Acceptable to Warning level
> 150	Unacceptable (incl. Curvilinear Data)



## Method Summary

SP1223	Sample 1		Sample 2		Sample 3		Sample 4	
Method Name	Mean	CV%	Mean	CV%	Mean	CV%	Mean	CV%
Abbott Enzymatic (n=2)	159.8	1.4	258.0	0.7	760.7	1.0	557.4	1.2
CLS Enzymatic (n=1)	162.6	0.0	262.6	0.0	747.3	0.0	549.7	0.0
Emit tox Salicylic Acid (n=9)	162.8	6.0	256.3	4.3	744.2	5.1	539.4	6.0
Roche Enzymatic (n=37)	160.1	2.8	252.5	2.5	704.6	3.3	537.1	2.4
Siemens Enzymatic (n=5)	150.6	4.3	244.3	4.1	699.7	1.6	513.1	1.9

*There are no Weqas or Participant supplied comments for Salicylate for this distribution.*

# Overall Performance

Traffic light colour	Performance category	Description
Red	Non-Return	(Any number of samples in distribution) and (No results submitted for any samples)
	Poor	(Any number of sample in distribution) and (one or more $ PI  > 3$ )
		(number of samples in distribution = 1) and ( $2 <  PI  \leq 3$ )
		(number of samples in distribution $> 1$ ) and (2 or more $ PI  > 2$ ) ( <i>Westgard <math>2_{2s}</math></i> )
		(number of samples in distribution $> 1$ ) AND (sign $x_i \neq$ sign $x_j$ ) AND $ x_i  > 2$ AND $ x_j  > 2$ * for at least one pair of scores $x_i$ and $x_j$ where $x_n$ with $i \neq j$ are the individual sample-analyte scores for the instrument in question ( <i>Westgard <math>R_{4s}</math></i> )
Orange	Warning	(number of samples in distribution $> 1$ ) and (single sample with $2 <  PI  \leq 3$ ) and (no sample with $ PI  > 3$ )
Yellow	Acceptable	(any number of samples in distribution) AND ( $1 < \text{Max }  PI  \leq 2$ )
Green	Good	$ PI  < 1$ for all samples Green is only assigned if all the individual samples are good.

Weqas  
Unit 6, Parc Tŷ Glas  
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750  
Email: [contact@weqas.com](mailto:contact@weqas.com)

# Weqas

GLOBAL PROVIDER OF QUALITY  
IN DIAGNOSTIC MEDICINE



EXTERNAL  
QUALITY  
ASSESSMENT



INTERNAL  
QUALITY  
CONTROL




REFERENCE  
MEASUREMENT  
SERVICES




EDUCATION &  
TRAINING

## Problem Solving and Cases

# Problem Solving Flow Chart






# Problem Solving Flow Chart

## IMPRECISION

# Problem Solving Flow Chart


START HERE






# Problem Solving Flow Chart

## INACCURACY



# Bias plot (1)

Sodium (mmol/l)



$$y = 0.9x + 9.6$$


$$r = 0.9995$$

$$IS = 5$$

$$Sy.x = 0.63$$

# Bias plot (1)

Sodium (mmol/l)




$y = 0.9x + 9.6$   
 $r = 0.9995$   
 $IS = 5$   
 $Sy.x = 0.63$

Imprecision – satisfactory  
Inaccuracy – identify error  
 $m = 0.9$ ,  $c = +9.6$  mmol/L – OK at 100 mmol/L  
2.5% negative bias at 130 mmol/L, 4% negative bias at 160 mmol/L  
Error – mixed. Two points calibration at 110 and 160 mmol/L  
**Cause – incorrect values for 160 mmol/L calibration**

# Bias plot (2)

Chloride (mmol/l)



$$y = 1.01x - 4.01$$


$$r = 0.9995$$

$$IS = 5$$

$$Sy.x = 0.49$$

# Bias plot (2)

Chloride (mmol/l)



$$y = 1.01x - 4.01$$

$$r = 0.9995$$

$$IS = 5$$

$$Sy.x = 0.49$$

Imprecision – satisfactory

Inaccuracy – identify error


$c = -4.0$  mmol/L

Error – systematic absolute. Results low by 4.0 mmol/L over whole range

**Cause – incorrect serum blank compensation**

# Bias plot (3)

Bicarbonate (mmol/l)



y = slope not calculated


$r = 0.8484$

IS = 1516

$S_{y.x} = 2.41$

# Bias plot (3)

Bicarbonate (mmol/l)




y = slope not calculated  
 $r = 0.8484$   
IS = 1516  
 $Sy.x = 2.41$

Imprecision – unsatisfactory,  $r = 0.8484$ ,  $Sy.x = 2.4$  mmol/L  
Not curvilinear  
Error – random  
**Cause – faulty syringe on instrument**

# Bias plot (4)

Calcium (mmol/l)



y = slope not calculated

r = 0.9810


IS = 190

Sy.x = 0.09



# Bias plot (4)

Calcium (mmol/l)




y = slope not calculated  
 $r = 0.9810$   
IS = 190  
 $Sy.x = 0.09$

Imprecision – unsatisfactory,  $r = 0.9810$ ,  $Sy.x = 0.09$  mmol/L  
Error – curvilinear data  
**Cause – time expired reagents**

# Bias plot (5)

ALP (IU/L)



$$y = 0.50x + 1.16$$


$$r = 0.9999$$

$$IS = 1$$

$$Sy.x = 1.30$$

# Bias plot (5)

ALP (IU/L)



$y = 0.50x + 1.16$   
 $r = 0.9999$   
 $IS = 1$   
 $Sy.x = 1.30$

Imprecision – satisfactory  
Inaccuracy – identify error  
 $m = 0.57$

Error – systematic proportional error. Results low by 41% over the whole range

**Cause – incorrect method group classification. The lab was using AMP not DEA buffer**

Weqas  
Unit 6, Parc Tŷ Glas  
Llanishen, Cardiff, CF14 5DU

Tel: 02920 314750  
Email: [contact@weqas.com](mailto:contact@weqas.com)



GLOBAL PROVIDER OF QUALITY  
IN DIAGNOSTIC MEDICINE



EXTERNAL  
QUALITY  
ASSESSMENT



INTERNAL  
QUALITY  
CONTROL



REFERENCE  
MEASUREMENT  
SERVICES



EDUCATION &  
TRAINING

## Case Study - example

### Summary Report

Analyte:

Overall Performance Category:

Running PI analyte –sample scores:

### Standard Report

Analyte results table:

Scoring Table:

Running PI scores:

Bias Chart (Absolute):


Bias Chart (Relative):

Precision Scores:

Sample Histograms:

Method Summary Data:

Programme: **Lipids** • Distribution Code: **LP0124**  
 Distribution Start: 22-Jan-2024 • Distribution End: 05-Feb-2024 • Report Issued: n/a • Report Status: n/a  
 Requested By: [gareth@weqas.com](mailto:gareth@weqas.com)



PI Scores	
Location	The Path Lab
Instrument Name	Cobas C 311
Instrument Serial #	2039-05
Cholesterol	Acceptable
Triglyceride	Acceptable
HDL Cholesterol	Poor
LDL Cholesterol	Good
<b>Overall % poor PI</b>	19%
<b>Overall % Non-return</b>	0%

PI Code	Meaning
N/A	Not enrolled for this analyte
NRR	Analyte enrolled but no results returned
N/S	This quantitative analyte not scored for this distribution

**\*\*\*\*\* Please note\*\*\*\*\***

The NQAAP expects 100% compliance on EQA returns. If you are enrolled for an analyte and have not returned a result (denoted by NRR in your PI table), this will be treated as a poor performance score. Please let us know as soon as possible if you are no longer providing this analyte as part of your diagnostic service.

**Distribution Specific Comment**  
 There is no specific comment for this distribution.

*No participant comments have been submitted for any instruments for this distribution.*

## Summary Report

Analyte: HDL

Overall Performance Category: Poor

Running PI analyte –sample scores: >50% good, some acceptable some poor. More poor and less good than median

**Programme:** Lipids • **Distribution Code:** LP0124  
 Distribution Start: 22-Jan-2024 • Distribution End: 05-Feb-2024 • Report Issued: • Report Status: n/a

**Participant Code:** WQ00942 • **Location:** The Path Lab • **Analyser Name:** Cobas C 311 • **Serial #:** 2039-05  
 Date Samples received: 23-Jan-2024 • Date of Analysis: 23-Jan-2024 • Operator Details: YW • Storage Conditions: Ambient Temperature (18 - 30 °C)  
**Analyte:** HDL Cholesterol • **Method:** Roche HDLC4 • **Kit:** HDLC4 (07528566 190)

HDL Cholesterol mmol/L						
	Your reported result (mmol/L)		Method: Roche HDLC4	Instrument Model: cobas c 311	Overall	Non-Scoring Reference Value
Sample 1	1.1	Mean	0.911	0.990	0.950	1
		SD	0.024	0.109	0.070	
		Uncertainty	0.0039	0.0789	0.0094	n/a
		n	60	3	86	
Sample 2	1.2	Mean	1.000	1.087	1.065	1.12
		SD	0.008	0.116	0.106	
		Uncertainty	0.0013	0.0840	0.0143	n/a
		n	60	3	86	
Sample 3	1.8	Mean	1.553	1.630	1.565	1.61
		SD	0.059	0.172	0.069	
		Uncertainty	0.0096	0.1238	0.0093	n/a
		n	60	3	86	
Sample 4	2.6	Mean	2.197	2.333	2.217	2.29
		SD	0.068	0.268	0.083	
		Uncertainty	0.0110	0.1934	0.0112	n/a
		n	60	3	86	



### Summary Report

Analyte: HDL

Overall Performance Category: Poor

Running PI analyte –sample scores: >50% good, some acceptable some poor. More poor and less good than median


### Standard Report

**Analyte results table:** your results higher than method mean, higher than instrument mean but n only 3. Non scoring ref value


### Scoring

Sample	1	2	3	4	Overall Performance
Target: (Method Mean)	0.911	1	1.553	2.197	
Weqas TAE	0.146	0.16	0.248	0.352	
PI	2.6	2.5	1.99	2.29	Poor


### Running PI



### Bias (Absolute)



### Bias (Relative %)



### Precision

	LP0823	LP0923	LP1023	LP1123	LP1223	This distribution: LP0124
Sy.x	0.02	0.05	0.04	0.05	0.06	0.03
IS	7	11	12	80	20	6

*Sy.x is the average deviation from the best fit line and is an index of scatter.*

IS Score	Interpretation
0 to 10	Good
11 to 150	Acceptable to Warning level
> 150	Unacceptable (incl. Curvilinear Data)

## Summary Report

Analyte: HDL

Overall Performance Category: Poor

Running PI analyte –sample scores: >50% good, some acceptable some poor. More poor and less good than median

## Standard Report

**Analyte results table:** your results higher than method mean, higher than instrument mean but n only 3. Non scoring ref value





**Scoring Table:** 3 samples poor PI, 1 sample acceptable but PI 1.99.

**Running PI scores:** LP0823 acceptable, since them got worse, mostly poor for LP0124

**Bias Chart (Absolute):** longstanding positive bias, this dist worse than previous, instrument mean close to method mean,  $y = 1.16 * x + 0.03$ , 16% positive proportional bias. Instrument SDs wide but  $n = 3$

**Bias Chart (Relative):** longstanding positive bias

**Precision Scores:** precision good, previously good or acceptable



■ All Results ■ My Method ■ My Instrument ● My Result

### Method Summary

LP0124	Sample 1		Sample 2		Sample 3		Sample 4	
Method Name	Mean	CV%	Mean	CV%	Mean	CV%	Mean	CV%
Elimination (n=8)	1.225	13.0	1.381	10.9	1.754	1.5	2.475	7.8
Enzymatic Immunoinhibition (n=14)	1.046	4.8	1.193	4.3	1.564	3.7	2.217	3.4
NMR Spectroscopy (n=1)	1.280	0.0	1.290	0.0	1.590	0.0	2.280	0.0
<b>Roche HDLC4 (n=60)</b>	<b>0.911</b>	<b>2.6</b>	<b>1.000</b>	<b>0.8</b>	<b>1.553</b>	<b>3.8</b>	<b>2.197</b>	<b>3.1</b>
Vitros (n=3)	1.033	6.3	1.200	0.0	1.633	4.0	2.333	2.8

There are no Weqas or Participant supplied comments for HDL Cholesterol for this distribution.

## Summary Report

Analyte: HDL

Overall Performance Category: Poor

Running PI analyte –sample scores: >50% good, some acceptable some poor. More poor and less good than median

## Standard Report

**Analyte results table:** your results higher than method mean, higher than instrument mean but n only 3. Non scoring ref value

**Scoring Table:** 3 samples poor PI, 1 sample acceptable but PI 1.99.

**Running PI scores:** LP0823 acceptable, since them got worse, mostly poor for LP0124

**Bias Chart (Absolute):** longstanding positive bias, this dist worse than previous, instrument mean close to method mean,  $y = 1.16 * x + 0.03$ , 16% positive proportional bias. Instrument SDs wide but  $n = 3$

**Bias Chart (Relative):** longstanding positive bias


**Precision Scores:** precision good, previously good or acceptable

**Sample Histograms:** samples 3 and 4 off the scale, samples 1 and 2 not gaussian distributions but our result to the right of majority

**Method Summary Data:** Roche method negative bias to other methods

# Problem Solving Flow Chart

## INACCURACY



# Problem Solving Flow Chart

## INACCURACY

