



Quantitative Faecal Hb EQA

Weqas
Unit 6, Parc Tŷ Glas,
Llanishen
Cardiff, UK
CF14 5DU
Tel: +44 (0) 2920 314750
E-mail: office@weqas.com

INTENDED USE

Weqas Quantitative Faecal Hb (Fhb) EQA samples are for in-vitro diagnostic use as an external quality assessment material for testing for the presence of Haemoglobin in faeces.

SUMMARY

External Quality Assessment (EQA), or proficiency Testing (PT) is an essential part of providing quality laboratory diagnostic services, and participation in EQA is required for laboratory accreditation to ISO 15189 and ISO 17025. EQA is the inter laboratory comparison and performance evaluation that extends throughout all phases of the healthcare diagnostic testing cycle.

PRODUCT DESCRIPTION

Material is dispensed aseptically into buffered vials specific to each instrument manufacturer. The base material is organic material which closely mirrors the basic constituents of human faeces. This material is then spiked with human whole blood (Hb) to cover the pathological and analytical range for Faecal Immunochemical tests (FIT).

STORAGE AND STABILITY

Samples are stable for 14 days in a refrigerator at 2-8°C.

The samples should be assayed on the day you receive them. If you are unable to assay the samples immediately upon receipt, please store at 2-8°C until analysis. Samples stored at 2-8°C prior to analysis must be brought to ambient temperature (18-30°C) and assayed within 30 minutes.

PROCEDURE

The samples should be treated the same as patient specimens and run in accordance with the instructions accompanying the test system being used.

Safely dispose of excess sample in accordance with local waste policy guidelines.

Always wear gloves to avoid contamination.

LIMITATIONS OF PROCEDURE

The Weqas Quantitative Faecal Hb (Fhb) EQA samples should be analysed according to the instructions within this document. If there is evidence of microbial or fungal contamination in the product, discard the vial.

The Weqas Quantitative Faecal Hb (Fhb) EQA samples require storage as described in STORAGE AND STABILITY and handling as described in PROCEDURE.

Accurate and reproducible results are dependent upon properly functioning instruments and reagents and the use of correct procedures.

ANALYTES COVERED

The Weqas Quantitative Faecal Hb (Fhb) EQA samples cover a relevant pathological and analytical range.

Analyte
Quantitative FIT

! CAUTION !

Human source material.

The base material is organic material spiked with human whole blood (Hb). Each unit of whole blood used in the manufacture of this product has been tested by CE marked IVD methods and found to be non-reactive for the presence of HBsAg and antibody to HIV-1/2, HCV and HIV-1 Ag. While these methods are highly accurate, they do not guarantee that all infected units will be detected. Because no known test method can offer complete assurance that the hepatitis B virus, hepatitis C virus, human immunodeficiency virus (HIV) or other infectious agents are absent, all products containing human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

Although every effort is made to ensure that the material is free from any known infectious agent, the samples should be handled as for clinical specimens.