


How to calculate your long term bias for your uncertainty calculation?

The following is an example on how to calculate the long term bias and between batch imprecision for serum sodium for samples analysed between 7/04/2015 and 14/01/2016. The information is available from your "End of batch report". These reports are listed at the end of your Distribution table and accessed via Lab Stats.



- Logout
- Home
- PDF Report Prefs
- Questionnaire

- ▶ Queued Jobs
- ▶ Distributions
- ▶ Users
- ▶ Pools
- ▶ Materials
- ▶ Analyses
- ▶ Regions
- ▶ Reports
- ▶ Methods
- ▶ Manufacturers
- ▶ Suppliers (C & R)
- ▶ Instruments
- ▶ Sections
- ▶ Scheme Admin
- ▶ Lab Orders
- ▶ Prices
- ▶ Audit
- ▶ Add Lab
- ▶ Add Region
- ▶ Mailer
- ▶ Utilities
- ▶ Report Statistics

Find Lab:

- ▶ Lab Stats
- ▶ Lab Profile
- ▶ Scheme Selection

List Distributions

Scheme: Mainline Chemistry

To view stats select a distribution from the list below.

Please note, during calculation runs for a distribution, reports are temporarily unavailable. If the distribution you are interested in does not appear on the list below, please try again later.

Code	Distribution Date	Return Date	Standard Report		Simplified Report	
			Summary Report	Full Report	Summary Report	Full Report
Mainline Chemistry						
QP	Mon 04/07/16	Mon 19/07/16	View on screen	View on screen	Request PDF	View on screen
OO	Mon 05/05/16	Tue 21/05/16	View on screen	View on screen	Request PDF	View on screen
QN	Tue 03/05/16	Tue 17/05/16	View on screen	View on screen	Request PDF	View on screen
QM	Mon 04/04/16	Mon 19/04/16	View on screen	View on screen	Request PDF	View on screen
QL	Mon 29/02/16	Mon 14/03/16	View on screen	View on screen	Request PDF	View on screen
CK	Mon 01/02/16	Mon 15/02/16	View on screen	View on screen	Request PDF	View on screen
QJ	Mon 04/01/16	Mon 19/01/16	View on screen	View on screen	Request PDF	View on screen
QI	Mon 30/11/15	Mon 14/12/15	View on screen	View on screen	Request PDF	View on screen
QH	Mon 02/11/15	Mon 16/11/15	View on screen	View on screen	Request PDF	View on screen
QG	Mon 05/10/15	Mon 19/10/15	View on screen	View on screen	Request PDF	View on screen

Page 1 of 23. [1](#) [2](#) [3](#) [4](#) [5](#) [6](#) [7](#) [8](#) [9](#) [10](#) [11](#) [12](#) [13](#) [14](#) [15](#) [16](#) [17](#) [18](#) [19](#) [20](#) [21](#) [22](#) [23](#)

End of Batch Reports						
Scheme	From	To	Created	Report	Request PDF	
Mainline Chemistry	NJ (06/07/2009)	NS (06/04/2010)	09/08/2011	Report	Request PDF	
Mainline Chemistry	NU (07/06/2010)	OC (07/02/2011)	15/08/2012	Report	Request PDF	
Mainline Chemistry	OE (04/04/2011)	ON (03/01/2012)	15/08/2012	Report	Request PDF	
Mainline Chemistry	OW (01/10/2012)	PF (01/07/2013)	17/03/2014	Report	Request PDF	
Mainline Chemistry	PQ (02/06/2014)	PZ (02/03/2015)	16/10/2015	Report	Request PDF	
Mainline Chemistry	QA (07/04/2015)	QJ (04/01/2016)	09/08/2016	Report	Request PDF	

Click on latest End of Batch Reports

Select the Scheme from the drop down menu


The following table will be provided for each measurand for each Section (analyser) registered with us.

Lab Code: AE . Section Name: Architect ci 16200. Scheme: Mainline Chemistry. Distribution Range: QA - QJ

Analyte: Sodium (mmol/l)

Method: Indirect ISE	M871	M872	M873	M874	M875	M876	M877	M878
Section Stats								
Mean reported results	106.4	112.8	120.1	127.5	134.8	141.8	149.0	155.0
SD reported results	0.9	0.8	1.2	0.5	0.5	0.8	1.0	0.7
CV(%) reported results	0.81	0.74	0.96	0.41	0.37	0.59	0.69	0.46
Number of results	5	5	3	4	4	5	5	5
Method Result Stats								
Mean method mean	107.7	114.2	120.8	127.4	134.1	140.7	147.5	154.2
Median CV	0.88	0.88	0.54	0.74	0.71	0.68	0.61	0.73
Overall Result Stats								
Median CV	0.92	0.88	0.59	0.75	0.71	0.71	0.66	0.73

CV Comparison




Lab Code: AE . Section Name: Architect ci 16200. Scheme: Mainline Chemistry. Distribution Range: QA - QJ

The End of Batch report

For a large number of our Schemes, material is prepared and distributed on more than one distribution, e.g. a batch of Serum chemistry samples consists of 8 levels and each level is distributed on 4 or 5 occasions over a 10 month period. This allows calculation of your overall bias and between batch imprecision over this period.

The report provides the mean, SD and coefficient of variation (CV%) of your results for each level over this time period. Reference Target values for your measurands (where provided) are available from your Standard reports.

Reference values for samples M874, M873, M877 and M 871 respectively



Sodium	Pools								
Method: Indirect ISE	M871	M872	M873	M874	M875	M876	M877	M878	Units
Your results									
Mean reported results	106.4	112.8	120.1	127.5	134.8	141.8	149.0	155.0	mmol/L
SD reported results	0.9	0.8	1.2	0.5	0.5	0.8	1.0	0.7	mmol/L
CV reported results	0.81	0.74	0.96	0.41	0.37	0.59	0.69	0.46	%
Number of results	5	5	3	4	4	5	5	5	
Reference Target	108.2	114	120	125.9	132.1	139.3	145.8	152.9	mmol/L
Method Result Stats									
Mean method mean	107.7	114.2	120.8	127.4	134.1	140.7	147.5	154.2	mmol/L
Median CV	0.88	0.88	0.54	0.74	0.71	0.68	0.61	0.73	%
Overall Result Stats									
Median CV	0.92	0.88	0.59	0.75	0.71	0.71	0.66	0.73	%


How to calculate bias

Compare your Mean reported results with the Reference targets or if these are not available your Mean method mean and calculate the linear regression equation of your results (y) on the Reference target (x). This will be available on the End of batch report on the next software release, however, meanwhile, this can easily be achieved in Excel.

Copy and paste the table into Excel and calculate $y=mx+c$. For the above example this produced the equation:

$$y = 1.1054x - 12.529$$

Mean reported results	106.4	112.8	120.1	127.5	134.8	141.8	149	155	mmol/L
Reference Target	108.2	114	120	125.9	132.1	139.3	145.8	152.9	mmol/L
Mean method mean	107.7	114.2	120.8	127.4	134.1	140.7	147.5	154.2	mmol/L



Decide on your critical points for expressing bias (x) e.g. sodium of 130mmol/L

Apply in the equation:

$$Y = 130 * 1.1054 - 12.529 = 131.17 \text{ mmol/L}$$

$$\text{Bias} = \frac{(x-y)}{x} * 100 \%$$

$$\text{Bias} = \frac{(130-131.17)}{130} * 100\% = + 0.9\%$$