

List of reference measurement services

This file was created on 02 November 2015 from the JCTLM-DB website (<http://www.bipm.org/jctlm/>)

Your search criteria: Reference measurement services; Analyte: creatinine; Analyte category: Metabolites and substrates; Matrix category: -

WEQAS, United Kingdom Phone : +44 (0)29 20 748355 Fax : +44 (0)29 20 748336 Web : http://www.weqas.com Contact person : Mr David Ducroq Email : dave@weqas.com	
Analyte	creatinine
Material or matrix	blood serum
Applicable material or matrix	serum
Quantity	Amount-of-substance concentration
Service measurement range	25 µmol/L to 600 µmol/L
Expanded uncertainty (level of confidence 95%)	1.24 %
Interlaboratory comparison results	RELA - IFCC External Quality assessment scheme for Reference Laboratories in Laboratory Medicine at http://www.dgkl-rfb.de:81/index.shtml
Measurement principle	ID-GCMS (Isotope Dilution Gas Chromatography Mass Spectrometry)
JCTLM reference measurement method/procedure	University of Ghent reference measurement procedure for creatinine