

Qualitative Simplified Participant Report Interpretation

Weqas
Unit 6, Parc Ty Glas
Llanishen
Cardiff
United Kingdom
CF14 5DU
Tel: 02920 314750
Email: weqas.enquiries@wales.nhs.uk
Programme Organiser: Mrs Annette Thomas

Participant Code
replaces Lab Code.

Participant Code: WQ00000

#####


Distribution details
and report status.

Programme: **phIGFBP-1** • Distribution Code: **PH0426**
Distribution Start: 22-Apr-2026 • Distribution End: 07-May-2026 • Report Issued: 15-Jun-2026 • Report Status: Re-Issued
Requested By: #####

Your % Poor PI (Performance Index) across all samples for all instruments, the Median All Participants % Poor PI and the 97.5th Centile % Poor PI.

This Distribution	
Your % Poor PI	0%
Median All Participant % Poor PI	0 %
97.5 th Centile % Poor PI	50%

Running PI analyte-sample Scores



PI Ranges	
0	Good
1	Acceptable
2	Warning
3, 4, 5	Poor

The running table shows the current distribution PI scores, and ultimately the previous 5 distribution PI scores for all analytes, for your Participant Code, plus all-Participant Median performance.

The PI Scores table shows performance for all analytes and all instruments.

PI Scores	
Location	#####
Instrument Name	#####
Instrument Serial #	
phIGFBP-1	Good
Overall % poor PI	0%
Overall % Non-return	0%

PI scores are converted into performance categories and colour coded.

Overall % Poor PI across all analytes for each instrument, and Overall % non return for each instrument.

PI Code	Meaning
N/A	Not enrolled for this analyte
NRR	Analyte enrolled but no results returned

Key for PI Codes.

***** **Please note** *****
Weqas are a UKAS ISO 17043 accredited EQA provider. Our accredited programmes are described on our UKAS schedule of accreditation which is available here: [WEQAS 17043 Schedule of Accreditation](#)

Distribution Specific Comment

There is no specific comment for this distribution.

Distribution Specific Comments from Weqas or participant.


No participant comments have been submitted for any instruments for this distribution.

Programme: phIGFBP-1 • **Distribution Code:** PH0426
Distribution Start: 22-Apr-2026 • Distribution End: 07-May-2026 • Report Issued: 15-Jun-2026 • Report Status: n/a

Participant Code: WQ00000 • **Location:** ##### • **Analysed Name:** ##### • **Serial #:** n/a
Date Samples received: 22-Apr-2026 • Date of Analysis: 01-May-2026 • Operator Details: n/a • Storage Conditions: 2 - 8 °C
Analyte: phIGFBP-1 • **Method:** Actim Partus • **Kit:** 31931ETAC / 31931ETAL (31931ETAC / 31931ETAL)

Scoring			
Sample	1	2	Overall performance
Your response	Positive	Positive	Good
Correct response	Positive	Positive	
PI	0	0	

Qualitative PI Ranges	
Score	Category
0	Good
1	Acceptable
2	Warning
3, 4, 5	Poor



Your overall performance for this analyte is:
Good

- There are no current performance alerts for this analyte

There are no Weqas or Participant supplied comments for phIGFBP-1 for this distribution.

Your reported response for each sample.

The correct response for each sample is shown here.

Your PI score is shown for each sample. Each PI score has been assigned by comparing your result with the correct response, taking into account sample concentration. Overall performance category is assigned using an algorithm based on the individual sample PIs.

Running Performance chart shows the PI (Performance Index) score for each sample for the current distribution and the last 5 distributions, with colour coding for each performance category.

The individual Location and instrument details will be shown here.

Overall Performance is stated here along with a basic interpretation of your performance and potential error identified.

Overall Performance is shown as a 'Traffic light', with performance colour coded to easily identify the performance category.

Analyte Specific Comments from Weqas or participant.