


Participant Code, Distribution details and report status.

Participant Code: **WQ00000**
 Programme: **POCT INR** • Distribution Code: **IN0922** • Analyte: **INR** • Units: **units**
 Distribution Start: 13-Sep-2022 • Distribution End: 11-Oct-2022 • Report Issued: n/a • Report Status: n/a

Analyte specific breakdown of performance categories for all samples.



	Overall		Sample 1	
	n	%	n	%
Good	12	28.6	12	28.6
Acceptable	6	14.3	6	14.3
Poor	5	11.9	5	11.9
Non-Return	19	45.2	19	45.2
Non-scored	0	0	0	0
Total	42		42	

Numerical analyte specific breakdown of performance categories for each sample and overall for your Participant, Group or Super Group dependent on access level and report chosen.

Method specific data table showing:


- Target Value to which the submitted result is compared against and Target Value type e.g. Method Mean.
- Mean and SD of your submitted results for this method / instrument.
- Method data (all participants within this method) is also supplied in the right-hand column.

If you have more than one meter type, this data will be shown in a separate table.



Method	Abbott I-Stat	
Target Value Type	Method Mean	
	Sample 1	
Target Value	2.704	
	My result(s)	Method
Mean	2.781	2.704
SD	0.456	0.29
Uncertainty	0.0996	0.0457
n	21	63

Histogram showing the spread of your results. Bars are colour coded to show performance categories across the range.



Histogram showing breakdown of results for each sample, showing the spread of results for your method, your instrument and your submitted results.

Weqas supplied comment for INR:
There is no comment for INR for this distribution

Distribution Specific Comments from Weqas or participant.